



246 Columbus Avenue  
 Roselle, NJ 07203  
 1.888.523.1600  
 1.908.259.9105 (Fax)  
 E-Mail: [comfortfitlabs@aol.com](mailto:comfortfitlabs@aol.com)

[www.comfortfitlabs.com](http://www.comfortfitlabs.com)

### DOCTOR & PATIENT INFORMATION

Doctors Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 City: \_\_\_\_\_ Account #: \_\_\_\_\_  
 Patient Name: \_\_\_\_\_  Male  Female Age: \_\_\_\_\_  
 Height: \_\_\_\_\_ Weight: \_\_\_\_\_ Shoe Size: \_\_\_\_\_  
 Shoe Type: \_\_\_\_\_ Shoes Enclosed:  Yes  No  
**Cast Enclosed for:**  LEFT  RIGHT  B/L  
**PLEASE MARK MEDIAL AND LATERAL MALLEOLI ON THE NEGATIVE CAST**

### CLINICAL INFORMATION

#### DIAGNOSIS:

**Stance Evaluation**  
 Calcaneus alignment to leg: \_\_\_\_\_° inverted or \_\_\_\_\_° everted  
 Leg alignment to floor: \_\_\_\_\_° varum or \_\_\_\_\_° valgum

### RICHIE BRACE® PRESCRIPTION

**RICHIE BRACE® (standard):** Full Flexion Ankle Hinge Pivot.  
 Can include enhancements for Posterior Tibial Tendon Dysfunction (check any or all):  
 Medial Heel Skive  4mm  6mm Navicular Accommodation  (please mark negative cast)  
 Adjust Limb Uprights for Tibial Varum  Yes  No (see measurements above)  
 FOR SEVERE PTTD, RECOMMENDED MEDIAL ARCH SUSPENDER (SEE BELOW)  
**SPECIAL MODIFIED VERSIONS OF STANDARD RICHIE BRACE®:**  
 **RICHIE SOCCER BRACE®** - Includes shin guard.  
 **LITTLE RICHIE BRACE®** - Pediatric application for shoe size 4 and under.  
 **RICHIE BRACE® RESTRICTED ANKLE PIVOT:** Limits ankle motion, yet allows smooth contact phase of gait.  
 Indications: DJD ankle & STJ, tarsal coalition, mild Charcot, lateral ankle instability, peroneal tendinopathy.  
**ENHANCEMENTS (optional):**  
 **MEDIAL ARCH SUSPENDER** – Adjustable lifting strap under talo-navicular joint for severe PTTD  
 **LATERAL ARCH SUSPENDER** – Adjustable lifting strap under calcaneal-cuboid joint for peroneal tendinopathy and severe lateral ankle instability.

**RICHIE BRACE® DYNAMIC ASSIST:** Full flexion pivot with spring hinges for dorsiflexion assist.  
 Patient requirements: 1. Dropfoot 2. Ankle dorsiflexion to at least 90° to leg 3. Stable knee **(must have all 3)**

**RICHIE BRACE® SOLID AFO:** Traditional full leg posterior shell w/balanced functional orthotic footplate.  
 Indications: Dropfoot with unstable knee, dropfoot with spasticity, Charcot Arthropathy.  
**STS Bermuda Casting Sock Required**

**RICHIE GAUNTLET®**  7"  9"  
 **RICHIE CALIFORNIA®**  
 GAUNTLET AND CALIFORNIA COLOR OPTION -  TAN  CHOCOLATE  
**Both The Richie Gauntlet and The Richie California require the STS mid leg sock**

### ALL RICHIE BRACES® HAVE THE FOLLOWING STANDARD FEATURES:

◆ Top Cover – Implus®	◆ Cover Length - Mets	◆ Limb Uprights Supports – Aligned Perpendicular to Foot Plate
◆ Color – Black	◆ Orthotic Foot Plate – Intrinsic Balance to Perpendicular	◆ Heel Stabilizer Bar - Included
◆ Heel Cup – 35mm		

COLOR OPTION -  FLESH TONE  WHITE

### RICHIE BRACE® MODIFICATIONS

**NOTE: NON-STANDARD BRACE MODIFICATIONS MAY HAVE EXTRA CHARGES – SEE PRICING SHEET**

Top Cover	Length	Heel Cup	Medial Heel Skive
<input type="checkbox"/> Implus (standard)	<input type="checkbox"/> to Mets (standard)	<input type="checkbox"/> 10 mm	For severe pronation control
<input type="checkbox"/> Spenco	<input type="checkbox"/> to Sulcus	<input type="checkbox"/> 14 mm	<input type="checkbox"/> 2mm
<input type="checkbox"/> EVA	<input type="checkbox"/> to Toes	<input type="checkbox"/> 18 mm	<input type="checkbox"/> 4mm
<input type="checkbox"/> Diabetic (Plastazote/Poron)	<input type="checkbox"/> add poron cushion to extension	<input type="checkbox"/> 35 mm (standard)	<input type="checkbox"/> 6mm

### CAST AND ORTHOTIC MODIFICATIONS

<input type="checkbox"/> Heel Lift _____ (inch)	<input type="checkbox"/> <b>Orthotic Plate Accommodation (please mark on cast)</b>	<b>Forefoot Posting</b> _____° Varus _____° Valgus
<input type="checkbox"/> Add Medial Arch Flange	<input type="checkbox"/> Navicular <input type="checkbox"/> Medial Fascia Band	<b>Note:</b>
<input type="checkbox"/> Add Lateral Clip	<input type="checkbox"/> Styloid 5 <sup>th</sup> Met <input type="checkbox"/> Other:	Not recommended as this will tilt entire brace to exact degree of posting.

#### SPECIAL INSTRUCTIONS:

**Accommodation location(s):**  
 (mark on illustration and on cast)

