



ComfortFit Orthotic Labs
NIGHT SPLINTS & WALKERS 2023

Products ideal for your patients who have Plantar Fasciitis, Achilles Tendonitis and other lower extremity overuse injuries.



Posterior Adjustable Night Splint *



Dorsal Night Splint

POSTERIOR ADJUSTABLE NIGHT SPLINT:

- Adjustable flexion straps provide gentle stretching of the plantar fascia and Achilles tendon
- Adjust from 10 degrees to 90 degrees dorsi-flexion
- Lightweight, sturdy, one piece flexible shell
- Cool foam liner, padded calf and ankle straps for additional comfort
- Shipped with toe wedge for additional stretch
- **Indications:** Plantar Fasciitis, Achilles tendonitis and other lower extremity overuse injuries.
- Available in 3 sizes (small, medium, large). *Extra Large available at \$19.50 /each

DORSAL NIGHT SPLINT:

- Light- weight, low profile design for increased comfort and compliance
- Slips on and adjusts easily
- **Indications:** Plantar fasciitis, Achilles tendonitis, drop foot, and post-static pain.
- Available in 2 sizes: small-medium and large-XLarge

L Code: L4396



Premium Air Walker

Available in low top or standard height

- Low profile design provides comfort & stability
- Air pump on front adjusts easily to achieve desired compression
- Adjustable medial/lateral air cell system inflates both sides equally to accommodate changes in edema during rehab
- Fixed, aluminum malleable uprights allow for the frame to contour better to the patient's leg
- Non-skid rocker sole is designed to help promote a natural gait
- **Indications:** Acute ankle sprains, soft tissue injuries, stable fractures, post-operative use.
- Available in 5 sizes

L codes: L4360



Premium Walker
Available in low top or standard height

- Low profile design provides comfort and stability
- Fixed aluminum malleable uprights allow for the frame to conform better to the patient's leg
- Non-skid rocker sole is designed to help promote natural gait
- **Indications:** Acute ankle sprains, soft tissue injuries, stable fractures, post-operative use
- Available in 5 sizes

L Code: 4386

ComfortFit Orthotic Labs
246 Columbus Avenue
Roselle, NJ 07203
1.888.523.1600
www.comfortfitlabs.com
contact@comfortfitlabs.com