



Prefabricated Orthotics for All Walks of Life

FUNCTIONAL ORTHOTICS FOR MEN'S AND WOMEN'S SPORT AND DRESS SHOES

COMFORT DRESS



A functional, intrinsically posted, full-width, 3/4-length, semi-rigid foot orthotic covered in a sulcus-length soft tissue top cover for easy fit and comfortable feel. This is a product of choice when a generic, cost-effective solution is desired in a shoe with limited room.

COMFORT ACTIVE



An orthotic for maximum control, extrinsically posted and covered in 1/8" of soft tissue top covering. With firmer shell and a deeper heel cup than our Dress model, this insole creates the perfect cost-effective sport device.

FOR DIABETIC OR SENSITIVE FEET

HEAT MOLDABLE



Heat-moldable insoles that provide maximum protection for the diabetic and sensitive foot. Full-length multi-density insole with a Plastazote top cover designed to be used with all diabetic shoe types. Medicare approved: suggested billing code A5512.

PEDIATRIC

COMFORT KID



Designed specifically for children. Multi-layer, firm anti-pronatory, firm medial and lateral flanges with soft foam top and anti-slip bottom. Children's 3/4 length. Sizes from Child's 3-4 to youth 7.

COMFORTFIT ORTHOTIC LABS

246 Columbus Avenue • Roselle, NJ 07203
Phone (888) 523-1600 • (908) 259-9100 • Fax (908) 259-9105
www.comfortfitlabs.com • contact@comfortfitlabs.com