



Sport-Fit Brace

L1902 Ankle Brace

The Sport-Fit Ankle Brace is our newest non-custom brace. The indications for this ankle brace are ankle sprains, recovery from ankle fractures, post-op, activity/sports-related support for chronic instability and prophylactic use. This is a commonly used brace for athletes to prevent ankle inversion/eversion. The brace can also be used on patients who have just come out of a short leg walker during their recovery process.



FEATURES

- 👉 Easy lace-up system
- 👉 Figure-eight nylon straps mimic athletic taping to help prevent inversion/eversion
- 👉 Elastic cuff at top of brace hides the laces and provides extra support
- 👉 Low-profile design fits into any athletic shoe
- 👉 Available in 6 sizes: XSmall, Small, Medium, Large, XLarge, XXLarge
- 👉 L Code: L1902

Date: _____ Account Number: _____

Account Name: _____

Address: _____

Item	Ankle Circ*	Item Number	Unit	Price	Quantity	Amount
XSmall	10" - 11"	US8832	each	\$		\$
Small	11" - 12"	US8833	each	\$		\$
Medium	12" - 13"	US8835	each	\$		\$
Large	13" - 14"	US8837	each	\$		\$
XLarge	14" - 15"	US8838	each	\$		\$
XXLarge	15" - 16"	US8839	each	\$		\$

*Measurement taken diagonally over the base of the heel and across the ankle bone

To order, please contact us by phone at 1-888-523-1600, by fax at 1-908-259-9105, or by email at contact@comfortfitlabs.com. You may view and download order forms and sell sheets on our website at www.comfortfitlabs.com.