

Torch Walker

Total Contact Orthotic Restraining Custom Hybrid

Features

- Triple density custom insert to relieve plantar pressures and increase comfort
- Soft foam padded collar decreases posterior proximal pressure
- Rocker sole to insure proper gait
- Integrated custom molded solid AFO

Indications

- Charcot deformities
- Chronic plantar ulcerations
- Non-union fractures
- Traumatic injury
- Partial foot amputations
- Equinus contractures
- Reflex sympathetic dystrophy
- Rheumatoid arthritis

Closures are laces, Velcro or any combination. You can order any height. The following are available: base and relief inserts, anterior shell tongue, amputation fillers, SACH heel, and integrated AFO.



To order please contact us by phone 1-888-523-1600 or by email at contact@comfortfitlabs.com. You may view download order forms and sell sheets on our website at www.comfortfitlabs.com

Codes

L-1960 A rigid orthosis to control the foot and ankle, custom molded from a model of the patient, custom fabricated, includes casting and cast preparation (if the device extends to a height no greater than 1.5 inches distal to the apex of the head of the fibula.)

Or

L-1940 A ankle foot orthosis, plastic or other material, custom fabricated (for device with a height below 1.5 inches distal to the apex of the head of the fibula.)

Plus the codes below

L-2330 A leather, or similar material, lacer molded from a model of the patient, that fits around the calf area, includes casting and cast preparation

L-2820 Addition to Lower Extremity Orthosis, Soft Interface for Molded Plastic, Below Knee Section

L-2232 Addition to lower extremity orthosis, rocker bottom for total contact ankle foot orthosis, for custom fabricated orthosis only

L-2280 Addition to lower extremity, molded inner boot

L-3020 Foot insert, removable, molded to patient model, longitudinal/metatarsal support each

L-5000 Partial foot, shoe insert, with longitudinal arch, toe filler (if the device includes an amputation filler code)

If patient is diabetic and qualifies can also use A5501 and A5503

NOT ALL ABOVE CODES APPLY TO ALL TORCH WALKERS ORDERED. BILL FOR WHAT YOU HAVE SUPPLIED. 1/23w